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EN | INSTRUCTIONS FOR USE

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ADITAMENTOS, COMPONENTES
E INSTRUMENTAL

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EN | INSTRUCTIONS FOR USE

Abutments, components and instruments

INSTRUCTIONS FOR USE

IMPORTANT INFORMATION

Read these instructions for use carefully before using the product.

This document contains basic information for the use of the prosthetic instruments and abutments or components of the NetPlus® original dental implant systems, from now on referred to as "NetPlus® products". This document has been written as a quick consultation guide for the practitioner responsible for the treatment, from now on referred to as "User", and does not provide sufficient technical specifications for the correct use of the NetPlus® products. It is not an alternative or substitute for specialised training and professional clinical experience.

The NetPlus® products must be used carrying out sufficient planning of the treatment and rigorously following the surgical and/or prosthetic protocols set by Biostar Laboratorios SLU. Read the specific surgical and prosthetic protocols for each product carefully, as well as the instructions for use and maintenance, before using any NetPlus® product. You can read them on our website, www.netplusdental.com, or request them from your official authorised Biostar Laboratorios SLU distributor.

1. Description of the system

The NetPlus® products are composed of dental implants, prosthetic components or abutments and surgical and prosthetic instruments.

Biostar Laboratorios SLU has created a wide range of dental implant systems for the carrying out of dentistry treatments in accordance with the existing scientific bibliography and current clinical standards. Biostar Laboratorios SLU develops, manufactures and commercialises latest-generation implant products and a wide range of restorative solutions to respond to the particular needs of each patient.

2. Description of the products

Instruments

The NetPlus® instruments are manufactured in stainless steel. Some instruments are laser-engraved to facilitate identification and use.

NetPlus® instruments

Insertion keys - Screwdriver tips.

Abutments

The NetPlus® prosthetic components or abutments are manufactured in grade 5 titanium, Cr-Co alloy and biocompatible plastic materials. Some NetPlus® components, abutments and screws have an anodised colour depending on the type of prosthetic platform to facilitate identification and use (see catalogue).

NetPlus® prosthetic components and components

Preparable straight abutments - Preparable angled abutments - Transepithelial abutments - Abutments for overdenture - Interphase abutments - Provisional abutments - Clinical screws - Cover screws - Healing abutments - Abutments with machined base - Impression abutments - Laboratory test abutments - Impression screws - Laboratory screws - Analogues - Scanbody - UCLA.

3. Indications for use

It is essential to previously evaluate the patient by means of a clinical and imaging diagnosis and study of models, as well as to determine the possible contraindications for the treatment, prior to carrying out a treatment with NetPlus® products. A correct treatment with dental implants must be based on a diagnosis, a treatment plan and a design of the prosthetic rehabilitation which meets the needs and expectations of the patient.

Instruments

The NetPlus® surgical and prosthetic instruments are designed as one of the fundamental parts of the implant treatment. With the instruments, you will be able to carry out all the necessary processes for the preparation of the bone tissue, the insertion of the NetPlus® dental implants and the manipulation of the different NetPlus® prosthetic components or abutments.

The NetPlus® instruments are re-usable.

Abutments

NetPlus® abutments and prosthetic components are designed to be connected to dental implants and are intended to support prosthetic rehabilitation. The prosthetic rehabilitation may be single or multiple, partial or total, and may be screwed or cemented to the abutments, or may be removable by means of implant anchorages or retainers.

The attachments or abutments in biocompatible plastic material (PEEK) are for provisional use and are designed to be used in the mouth for a period of time no longer than 180 days. The removable prostheses with anchorages or retainers must be monitored at least once per year by the clinical professional. The correct intraoral adjustment of the prosthesis and the retention between the anchorages or retainers and the dental prosthesis must be verified. A loss of adjustment and/or a poor retention of the prosthesis put the prosthetic rehabilitation and the dental implants in danger.

The NetPlus® abutments are single-use products.

4. Intended patients

Abutments

The target population for these medical devices are those patients who, in the opinion of the clinician assessing their case, require treatment in the mandible and maxilla, and are intended to restore chewing function, as well as providing an improvement in occlusion and the distribution of occlusal loads, phonation, aesthetics, prosthetic stability and maintenance of the alveolar process, as well as a positive psychological effect.

These patients must have completed their growth. The age guideline is older than 18 years, this can be confirmed after cessation of growth with a wrist X-ray. Patients at an early age patients will be specifically informed of the changes they may need to make to the prostheses during their lifetime with a with the consequent aesthetic and economic impact due to functional ankylosis of the implant in living tissue.

Instrumental

The target population for these medical devices are those patients who, in the opinion of the clinician assessing their case, require treatment in the mandible and maxilla, and are intended to restore chewing function, as well as providing an improvement in occlusion and the distribution of occlusal loads, phonation, aesthetics, prosthetic stability and maintenance of the alveolar process, as well as a positive psychological effect.

The target population age, approximately, can be from 6 years old with orthodontic treatment.

5. Contraindications

It is necessary to carry out a preoperative medical examination of the patient and to determine if absolute and/or relative contraindications or risk factors for the treatment with dental implants exist. The user is responsible for evaluating the potential risks and benefits of the implant treatment and for taking the pertinent decisions depending on the case. The NetPlus® products must not be used on patients who lack the minimum medical conditions for receiving a prosthetic treatment and rehabilitation with dental implants.

Relative contraindications

Age of the patient - Uncontrolled systemic disease - Poor oral hygiene - Active local infection - Pregnancy - Therapy with bisphosphonates - Sjogren syndrome - Uncontrolled endocrine disease - Radiation to head and neck - Allergy to stainless steel. - Stress - Psychological factors - Valve prostheses - Uncontrolled periodontal disease - Smoking - Alcoholism - Drug addiction - Addictions in general.

Absolute contraindications

The manufacturer does NOT recommend the implant treatment in the following cases: Cardiovascular diseases or disorders - Haemophilia - Cirrhosis of the liver - Terminal pathologies - Epilepsy - Neoplasms under care with chemotherapy - Blood dyscrasias - Cancer - Uncontrolled diabetes - Untreated periodontitis - Infectious diseases - Leukaemia - Alterations of the immune system - Treatment with immunosuppressives - Allergy to titanium - Active treatments with radiotherapy - Indefinite treatments with corticosteroids - Treatments with anti-coagulants - Involuntary muscle disorders - Serious liver disorders - Endocrine disorders - Psychiatric and/or psychological disorders.

Biostar Laboratorios SLU will update these contraindications when necessary, based on the information obtained through clinical studies and post-sales information.

6. Storage

The NetPlus® products must be stored in their original packaging, in a clean, dry place, protected from sunlight and adverse conditions.

7. Packaging and supply state

The NetPlus® instruments and prosthetic components or abutments are supplied packaged in a thermosealed plastic bag; the product inside is identified with a label stuck to the back.

Included in the abutments, the removable labels serve as a traceability control of the abutment. This is achieved by sticking one label to the implant card for the patient, and the rest for the clinic's record.

Biostar Laboratorios SLU recommends that the traceability of the abutments used be maintained by recording the reference and the lot used in in the patient's file.

Biostar Laboratorios SLU guarantees that all its products follow an exhaustive process of manufacture, control and cleaning before being packaged. The correct integrity of the packaging must be controlled before using the NetPlus® products, checking that it is not damaged. If there is any defect or damage, it must not be used and must be notified to Biostar Laboratorios SLU or to the official authorised distributor.

8. Cleaning and sterilisation

The instruments and prosthetic components or abutments of the NetPlus® dental implants systems are supplied NON- sterile.

The instruments and the prosthetic components or abutments which are to come into contact with the mouth must be cleaned, disinfected and sterilised before use by following the cleaning, disinfection and sterilisation protocol for the Biostar Laboratorios SLU instruments, surgical boxes and abutments. You can consult it on our website www.netplusdental.com or in the NetPlus® catalogues.

9. Precautions

Surgery and Instruments

The surgical protocols describe in detail the precautions and important factors to take into account during the surgical process of fitting the dental implants. Prior to the intervention, it is essential to have consulted the surgical protocol of the product, as well as all the indications about its correct usage indicated in the catalogues of the NetPlus® products.

Particular attention must be paid and the necessary measures taken to prevent the instruments, of such small size, from being swallowed or breathed into an airway. All the NetPlus® instruments for manual use incorporate retaining elements between each other or for use with dental floss, to avoid accidental ingestion.

It is important to reduce the time of surgery and the damage done to the area receiving the implants, paying special attention to bone trauma and the asepsia of the surgical field. The correct preparation of the surgical site is a fundamental factor in the success of the treatment. This process requires that the surgical protocol for drilling set by the manufacturer for each type of implant and the correct use of the surgical drills, both in terms of the number of uses and of the working speeds, be followed rigorously. The surgical drills and taps must be used with constant abundant irrigation to prevent overheating of the bone.

Take special care not to apply excessive forces or torques during the insertion of the implant. It is important to control the insertion torque of the implant, both in insertion with contra- angle, and with insertion with ratchet wrench.

Attention

When a surgical instrument with a limited number of uses is used, that number must not be passed, as this would put the process of osseointegration, and therefore the treatment, in grave danger.

The maximum torque for insertion of the dental implants is 50 Ncm. Exceeding the maximum insertion torque indicated for the implants can cause serious damage to the dental implant, to its connection, to the transporter and to the clinical screw included.

Prosthetic rehabilitation

The surgical protocols describe in detail the precautions and important factors to take into account during the surgical process of fitting the prosthetic rehabilitation. The design, the type of rehabilitation and the execution of the dental prosthesis must be studied prior to the insertion of the dental implants. Take special care not to apply excessive forces during the tests of the abutments or dental prostheses; this must always be done manually. Pay special attention to the definitive adjustment of the rehabilitation; the torque applied to the definitive clinical screws must be correctly controlled following the manufacturer's indications, whether carried out mechanically or manually.

NetPlus® Torque table for abutments

Refer to the recommended torques for each system and abutments in the table of torques on the website, www.netplusdental.com.

Some instruments are designed for use with a contra-angle. In order to guarantee a correct connection of both products, please see the contra- angle manufacturer's instructions for use.

Attention

Applying insufficient torque can cause incidents of loosening of said screws. If this occurs, we recommend replacing them with new products, as retightening a screw increases the risk of fatigue and fracture of it.

Exceeding the maximum torque applied to the clinical screws indicated by the manufacturer can cause serious damage to the dental implant, to its connection, to the prosthetic abutment and to the clinical screw. Always check the passive adjustment of the implant- prosthesis unit. An incorrect passive adjustment puts the prosthetic rehabilitation in grave danger and can cause future damage to the dental implant, to its connection, to the prosthetic abutment, to the clinical screw and/or cause the loss or reabsorption of bone.

10. Warnings

The planning of the treatment and the fitting of the NetPlus® dental implants requires specific dentistry training. The specifications of the product alone do not guarantee a proper use of it. The users are recommended to attend specialised theoretical/practical training courses to become familiar with the work techniques and protocols suitable for each product, including biomechanical, radiographical and prosthetic requirements associated with the implant treatment. Biostar Laboratorios SLU periodically holds diverse courses so that the users can update their knowledge of the work procedures and techniques of the NetPlus® implants systems. It is the user's obligation to inform themselves of and familiarise themselves with the current state of technique for each product and of its possible applications.

Before using NetPlus® products, it is necessary to know and be familiar with the corresponding surgical and prosthetic procedures. The patient must meet sufficient anatomical and psychological conditions to undergo the treatment with dental implants. The user must ensure the correct planning of the treatment, guaranteeing its safety with a minimum margin for error and respecting the vital oral structures and the general health of the patient.

The procedure for use of the surgical drills, taps and other surgical instruments needed for the fitting of the implants is specified in the surgical protocols corresponding to each type of implant. The fitting of the implant and the prosthetic planning must adapt to the patient's individual conditions, and in particular to a correct distribution of forces. The passive adjustment in the prosthetic rehabilitation structure and the adjustment of the occlusion must be achieved, and the appearance of excessive lateral forces avoided. A treatment with an insufficient number of implants, or an unsuitable choice of size and position to withstand and transmit the loads expected, can cause mechanical failure of the implant, of the abutment or prosthetic components or of the clinical screws due to overburdening or fatigue and substantial bone loss around the implant.

A qualitatively or quantitatively insufficient alveolar bone, the appearance of infections or diseases in general, and changes in the habits of the patient, are some potential causes of failure of the osseointegration and of the treatment. The lack of bone tissue or soft tissue can give rise to a deficient aesthetic result. An unsuitable prosthetic rehabilitation can cause the failure of the whole rehabilitation treatment.

The reuse of single- use products carries the danger of possible deterioration of their properties. The geometry of the product may have suffered damage during previous use or due to an unsuitable use. This implies the risk of surgical or prosthetic failure and/or harm to the patient's health.

Electrosurgery is not indicated due to the conductivity of the dental implants.

11. Adverse or side effects

Many adverse or undesired effects can present themselves in current implantology. These are documented in the specialised scientific bibliography and published in journals and books specialised in the dentistry sector.

However, the most relevant are:

Post-operative discomfort - Localised pain - Localised inflammation - Haematomas - Haemorrhages - Localised or systemic infections - Difficulty speaking - Bone loss and/or fracture - Loss of the implant - Damage to adjacent teeth - Fractures to the implants and/or prosthetic components- Fracture of instruments- Injury to adjacent anatomical structures (inferior alveolar nerve or maxillary sinus).

Biostar Laboratorios SLU will update these side effects when necessary, based on the information obtained through clinical studies and post- sales information.

12. Information on responsibility, safety and guarantee

The indications for use and handling of the NetPlus® products are based on the international published bibliography, the current clinical standards and our clinical experience with our products, and therefore should be understood as indicative general information. The handling and use of the NetPlus® products, being out of the control of Biostar Laboratorios SLU, are sole responsibility of the user of them. Biostar Laboratorios SLU, its subsidiaries and/or official authorised distributors assume no responsibility, whether express or implicit, full or partial, for the possible harm or damages which may result from the improper handling of the product or from any other event not specified in its protocols and manuals for the correct use of its products.

The user of the product must ensure that the NetPlus® product used is suitable for the intended procedure and purpose. Neither these instructions for use, nor the work and product handling protocols, absolve the user of this responsibility. The use, handling and clinical application of the NetPlus® products must be carried out by qualified professional staff, with the necessary qualifications in compliance with the legislation in force in each country. The use, handling and/or application of the NetPlus® products, whether fully or in part, at any stage, by unqualified staff or staff without the necessary qualifications, automatically renders any kind of guarantee void and can cause serious harm to the health of the patient.

The NetPlus® products form part of their own system, with their own design characteristics and work procedures. This includes the dental implants, abutments

or prosthesis components and the surgical or prosthetic instruments. The use of NetPlus® products in combination with elements or components from other manufacturers can cause failure of the treatment, cause tissue damage, cause damage to the bone structures, inadequate aesthetic results and serious harm to the patient's health. For this reason, only original NetPlus® products must be used.

The clinical professional responsible for the treatment is the sole person responsible for ensuring the use of original NetPlus® products and for using them in conformity with the corresponding instructions for use and handling protocols for the whole duration of the implant treatment process. The use of components, instruments or any other non- original NetPlus® product used alone or in combination with any of the original NetPlus® products will automatically void any kind of guarantee of the original NetPlus® products. Consult the Biostar Laboratorios SLU Guarantee Programme.

Any problem or incident with regard to the device should be reported to the manufacturer, to which the user should also send the affected product. In the event of a serious incident, the user should also submit a report to the manufacturer, as well as the relevant competent authority, according to the requirements of local legislation

The SSCP is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI.

Warning

Not all the NetPlus® products are available in all countries. Consult their availability.

The NetPlus® brand is a registered trademark, as are some of the products mentioned or not mentioned in this document. Biostar Laboratorios SLU reserves the right to modify, change, remove or develop any of the products, prices or technical specifications referenced in this document or in any of its catalogues without prior notice. All rights reserved. The whole or partial reproduction or publication of this documentation, by any means and in any format, is forbidden without the corresponding authorisation in writing from Biostar Laboratorios SLU.

13. Magnetic resonance

The safety and compatibility of the NetPlus® products in a magnetic resonance environment has not been evaluated, nor have they undergone testing to see if they heat up or migrate in said environment.

14. Removal of the material

The dental implants, abutments and kit, as well as each product and consumable used during surgery for placing dental implants and which are disposed of may pose a health hazard to people handling them. For this reason, you are advised to check and comply with current legislation on biohazardous waste.

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EN | CLEANING, DISINFECTION AND STERILISATION

Cleaning, disinfection and sterilisation

The protocols described in this section must only be carried out by personnel qualified to clean, disinfect and sterilise the dental materials specified here in.

CLEANING AND DISINFECTION INSTRUCTIONS

Applicable for instruments, surgical and prosthetic boxes.

Disassembly

1. Dismount* the appropriate instruments, for example manual ratchets, drills or drill stops.
2. Remove the various components from the surgical or prosthetic box for correct cleaning.

Cleaning and disinfection

For disinfecting instruments and surgical boxes:

1. Submerge the instruments in a detergent/disinfectant solution** suitable for dental instruments to help eliminate any adhered biological residues. If an ultrasound bath is available***, confirm that the detergent/disinfectant solution is indicated for use with this type of equipment.
2. Manually remove any biological residues with a non-metallic brush and pH-neutral detergent.
3. Rinse with copious water.
4. When cleaning the surgical and prosthetic boxes, always use a pH-neutral detergent and non-abrasive utensils to avoid damaging the surface of the boxes.
5. Dry the materials with disposable cellulose, lint-free clothes or compressed air.

For disinfecting plastic caps and spacers:

1. Submerge in a neat benzalkonium chloride solution for 10 minutes.
2. Rinse with distilled water.
3. Dry the caps and spacer before use.

Inspection

1. Check that the instruments are perfectly clean; if not, repeat the cleaning and disinfection steps.
 2. Discard any instruments with imperfections and replace them before the next procedure.
 3. Check that the instruments and the surgical and prosthetic boxes are perfectly dry before reassembling the parts and proceeding to their sterilisation.
- * See assembly and disassembly instructions.
- ** Follow the instructions from the disinfectant's manufacturer to determine the correct concentrations and times.
- *** Follow the instructions from the ultrasound bath's manufacturer to determine the correct temperature, concentration and times.

STERILISATION INSTRUCTIONS FOR STEAM AUTOCLAVE

Applicable to orthodontic implants, abutments, instruments, surgical and prosthetic boxes, pins, fixation screws and laminates.

1. Introduce each material separately in individual sterilisation bags, then seal the bags. For joint sterilisation, place the instruments in their surgical box, introduce the box into a sterilisation bag and seal the bag.
2. Place the bags to be sterilised in the autoclave.
3. Sterilise in a steam autoclave at 134°C/273°F (max. 137°C/276°F) for 4 min (minimum) and at 2 atm. Torque wrenches must be sterilised in 3 vacuum cycles at 132°C/270°F for a minimum of ≥ 4 minutes and vacuum-dried for a minimum of 20 minutes.

For the United States only: The validated and recommended sterilisation cycle for the US must be performed in a steam autoclave at 132°C/270°F for at least 15 min and with the drying time of at least 15- 30 min.

Important:

- Make sure the drying stage is allowed to run to completion, otherwise the products may be damp.
- Check the sterilisation equipment if the materials or sterilisation bags are damp at the end of the sterilisation cycle.
- Perform the necessary maintenance actions on the autoclave according to the established periodicity and following the manufacturer's instructions.

STORAGE OF NETPLUS® PRODUCTS

- Store the products in their original packaging and in a clean, dry location until they are used.
- After sterilisation, keep the products in the sealed sterilisation bags and in a clean, dry location.
- Never exceed the use by date indicated by the manufacturer of the sterilisation bags.
- Always follow the indications of the manufacturer of the sterilisation bags.

GENERAL RECOMMENDATIONS

- Never use damaged or dirty material; never reuse single-use products. The user is responsible for following the instructions described in this document correctly.
- The attention to piercing or sharp elements. Gloves should be worn when cleaning the materials to avoid accidents during handling.
- Follow the safety instructions indicated by the manufacturer of the disinfectant agent.
- The product's sterility cannot be guaranteed if the sterilisation bag is open, damaged or damp.
- Respect all stages of the sterilisation process. If the materials or sterilisation bags contain traces of water or moisture, check the autoclave and repeat the sterilisation.
- Orthodontic abutments and implants are supplied UNSTERILISED and must always be sterilised before use.
- Instruments and surgical and prosthetic boxes are supplied UNSTERILISED and must always be sterilised before use and cleaned and disinfected after use.
- The sterilisation, cleaning and disinfection processes gradually deteriorate the instruments. Inspect the instruments thoroughly to detect any signs of deterioration.
- Avoid contact between products made from different materials (steel, titanium, etc.) during the cleaning, disinfection and sterilisation processes.
- Biostar Laboratorios SLU recommends these instructions are implemented for the correct maintenance and safety of their products; accordingly, the company refuses any liability for any damage to the products that could arise if the user applies alternative cleaning, disinfection and sterilisation procedures.

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